


PH-28 | 41535-000

SPECT QA Phantom JSP



For daily quality control in SPECT and PET imaging




Optional parts: holder

SHOW MORE!



outer phantom



line source phantom



cold spot phantom



hot spot phantom



geometric distortion phantom



dose linearity phantom

FEATURES

A set of test units for daily QA of SPECT/PET

APPLICATIONS

SPECT and PET

EVALUATION PARAMETERS

- | Uniformity
- | Spatial resolution
- | Dose linearity
- | Image distortion

DESCRIPTIONS

SET INCLUDES

- | | |
|--------------------------------|------------------------|
| 1 outer phantom | 1 petroleum jelly |
| 1 line source phantom | 1 screwdriver |
| 1 cold spot phantom | 3 kinds of extra screw |
| 1 hot spot phantom | 1 Injection needle |
| 1 dose linearity phantom | 1 storage case |
| 1 geometric distortion phantom | manual |
| 1 phantom holder | |

MATERIALS

Phantom: methacrylic resin

SPECIFICATIONS

Phantom size:
22dia. x 22 cm/ 8.7dia. x 8.2 in

COMPILES WITH

JIS Z 4922



41535-010

Holder and accessories

Optional Parts for PH-28 and 30

Specify the manufacturer and type of the scanner

